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**INTENDED USE:**

The product is used for the transportation of patients in emergency or in situations demanding a quick and safe evacuation of the patient. Specially designed for confined spaces such as narrow corridors or revolving stairs with limited access. It has to be used exclusively by professional and qualified personnel, duly trained for emergency care. Can be used in ambulances and in hospital and sports environments. The stretcher can be converted from a transport chair to a wheeled stretcher for patient transport and in horizontal position.

**CAUTION (!):**

- In case of damage or malfunction of your product DOT NOT use the product and please contact the authorized and trained representative.
- Do not use the product for a purpose not indicated in this manual.
- ΜΟΒΙΑΚ S.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration on the frame of the product.
- ΜΟΒΙΑΚ S.A reserves the right to change the information contained in this document without previous notice
- There is a  $\pm 3\%$ . Deviation for the Dimensions of the Products.
- It is always recommended to follow the procedures designed for emergency response.
- Improper use of the device can cause injury to the patient.

**GENERAL WARNINGS:**

- If you don't read this user manual, it is preferable not to use this device or another available part.
- If you don't understand the usage or the precautions, please contact with the dealer or the proper technical person before you use the device because may be caused damage.
- READ carefully the entire manual before using the product.
- Save this document for future reference.
- The product should be used only by qualified personnel.
- Keep away the product from heat sources.
- Be careful when children are nearby and do not allow the children to play or operate with the device.
- Special care is always required where there are moving parts that can cause limb entrapment or injury.

**BEFORE EVERY USE:**

1. Check if the frame is damaged in order to guarantee a safe use of the product (there are no cracks or fractures in the frame).
2. In case of damage do not use the product and contact your dealer for further instructions.
3. Check, if the product is properly assembled and all screws are secure and well screwed on.
4. Always check the condition of wear for the mechanical parts to ensure the absolute safe use of the product for people and objects.
5. Periodically check its sturdiness and the stability.
6. Check for any tears in the surface fabric of the stretcher.

**USE OF PRODUCT:**

Open the stretcher.

Successfully move the patient on the stretcher and secure with the safety belts.

**CLEANING AND DISINFECTION:**

- Use only a wet cloth and mild soap on surfaces where they need cleaning. Then dry thoroughly before you use the product.
- Dry away from direct sunlight.
- DO NOT leave the product wet.
- Do not use chemical cleaners to clean the frame. This could cause damage to the surface of the product.
- Never use acids, alkalis or solvents such as acetone
- If you need to disinfect the product, use a common, mild disinfectant.
- Protect the item from scratches, cuts and punctures.
- Can be sterilized.

**DECLARATION OF CONFORMITY:**

We are solely responsible for declaring that the Medical Devices mentioned in this statement are of Low-Risk Class (Class I) and comply with the requirements of the European Regulation 745/2017 and where appropriate, the standards and legislation referred to.

**ENVIRONMENTAL PROTECTION:**

If one day you find that your product needs to be replaced or it is no longer working for you, consider protecting the environment:

- 1) Do not dispose your product along with the rest of the public waste (this is also the meaning of the shown recycling sign).
- 2) Contact your Public Authorities and they will instruct you of the Recycling centers to which your product must be disposed.
- 3) Correct disposal of your product helps the protection of the environment as well as the recycling of the product's components.

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## WARRANTY:

The product is guaranteed for a period of two (2) years from the purchasing date. The warranty apply on factory defects & DOES NOT apply on damages resulting from misuse, poor maintenance, modification, overuse, or non-compliance with the product's operating instructions. Also DOES NOT apply on the parts, that can wear out during use or over time. Damages or defects caused by natural disasters, unauthorized maintenance or repair, power supply problems (where applicable) & carriers are not covered by the warranty. No service or spare part is applied unless it has been inspected by an authorized service department or the manufacturer of the product. Shipping costs for goods of warranty are to be paid by the end user- customer or by dealer.

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DIMENSIONS HORIZONTAL POSITION	DIMENSIONS SITTING POSITION	DIMENSIONS CLOSED	CARTON DIMENSIONS
192 x 50 x 20 cm	BACKREST HEIGHT: 137 cm SEAT HEIGHT (FROM FLOOR): 52cm SEAT & BACKREST WIDTH: 50 cm	98 x 54 x 18 cm	98 x 54 x 18 cm

MBK CODE	PIECES PER CARTON	PIECE G.W.	PIECE N.W.	CARTON G.W.	CARTON N.W.
0805919	1	11 KG	9.6 KG	11 KG	9.6 KG

- Lightweight, easy to carry.
- Safe to use and easy to sterilize.
- Includes two orange safety belts. (code 0805927)
- With two wheels of diameter 150mm
- Foldable
- Made of high strength aluminum and waterproof PVC fabric.
- Can be converted from a transport chair to a wheeled stretcher for patient transport and in horizontal position.

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## INCIDENT REPORTING:

Any serious incident that occurs in relation to use of the device must be reported to the manufacturer and a competent authority of the Member State in which the user and / or patient resides.

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USER MANUAL: TRANSPORT STRETCHER CHAIR (2 IN 1)

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MOBIAK S.A