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Technical aid to support mobility. It serves to facilitate walking, increasing stability with a greater support base and a safer movement. The crutches are indicated to be used by people who have functional deficits in one or both lower extremities caused by pathologies such as arthritis and osteoarthritis, neurological diseases, surgeries or traumatic injuries. The use of the crutches allows the user to maintain independence to move around and helps to have a good, dynamic and active quality of life.

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If the crutches are not adjusted to the correct height, it can create bad posture vices that over time can cause injuries and muscle contractures in the back, arms and hands.

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Our products are guaranteed against any manufacturing defects for two years from the date of the invoice. For more information please refer to our general sale conditions.

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Any serious incident related to the product must be reported to Sucosores de Manuel García, S.L. and to the competent authority of the Member State in which the user and / or patient are established.

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1. The crutch height should be adjusted to suit each user's height in order to guarantee their safety and comfort.
2. To know the perfect height of the crutch you need to make sure that the forearm collar has to be 5cm below the height of your elbow.
3. To adjust the height, remove the fastener to relieve the mobile extension of the crutches.
4. When desired height is reached, lock tube into place by aligning both holes and introducing the fastener.
5. To be sure tube is locked into place by pushing down on or by pulling tube.

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1. The very first thing you should do is to select the most convenient and suitable crutches for you.
2. The crutch must be held in the opposite hand to the affected leg, either from weakness or surgery or you can use two crutches.
3. It is very important to check the height of the crutch, since it has to be adjusted according to the user's height.
4. The elbow should be slightly bended when you are up supported by the crutch.
5. The handle of the crutch must be at the same level as the wrist and the greater trochanter (hip bone).
6. There are different types of crutches handles. You should choose the one that provides you the most comfort.

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Please follow your therapist's advice before following the below advice:

- The posture when standing should be completely straight with a firm grip on the crutch.
- When stepping forward with your weaker leg, swing simultaneously the crutch in front of you in the same distance, so that the tip of the crutch and your front foot are even.
- Press down the crutch to remove some of the pressure on your weak leg.
- When you step with your strong leg, be sure to reach beyond the crutch.
- To turn, do it on your strong leg.

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When using two crutches we can use the last method of walking alternating both legs and both crutches.

If it is a major injury we can use the following method;

1. The first move it will be a small step with both crutches keeping the distance in between.
2. Then step forward with the injured leg until foot is in line with the crutches.
3. Press down the crutches to relieve pressure on the injured leg and to keep the balance with the crutches.
4. Lastly step forward with the strong leg to reach the other leg.

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1. Hold the crutches or the rail if you are using just one crutch.
2. Use the stronger leg to perform the first step
3. Next, bring the injured leg and the crutches to the same step.

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1. Put your feet near the steps and place your crutch on the first step down.
2. Step down first with your weak leg. Bring your stronger leg to the same step.