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PURPOSE OF USE:

This product allows third parties to wash the hair of people who are bedridden without wetting him/her or the bedding sheets. It can also be used by anyone and as a portable bath.

CONTENTS:

1. Inflatable basin with hose for easy water drainage.
2. Plastic tube & cord.

Water Basin

Water Bag, optional

CHARACTERISTICS:

1. Bathroom basin, double reinforced for greater durability & comfort. Also available with drain valve & hose for easy & quick water release. With water bag in case of inaccessible access to water source, **optional (0806232) with capacity of 8L.**

PREPARATION INSTRUCTIONS FOR USE:

1. When inflate the basin, close the valve immediately since it is not returnable valve. When deflate the basin, pinch the air nozzle with your fingers to deflate air.
2. Place the basin on a flat surface.
3. Place the drain hose in a basin or sink, make sure it is lower than the basin to prevent water from flowing.

4. Open the drain valve while using the basin to drain the water.

5. Place the basin under the user's neck adjust until he or she feels comfortable.

NOTE:

The water bag illustrated is **optional** (0806232)

Water bag, optional

ENVIRONMENTAL PROTECTION:

If one day you find that your product needs to be replaced or it is no longer working for you, consider protecting the environment:

- 1) Do not dispose your product along with the rest of the public waste (this is also the meaning of the shown recycling sign).
- 2) Contact your Public Authorities and they will instruct you of the Recycling centers to which your product must be disposed.
- 3) Correct disposal of your product helps the protection of the environment as well as the recycling of the product's components.

CLEANING INSTRUCTIONS AND WASHING INSTRUCTIONS

- Wash the basin with mild, antiseptic soap.
- Clean it with a wet cloth, rinse with clean water, and then dry with a dry cloth.
- Do not use detergents or other solvents.
- Do not wash it in the washing machine and do not use a dryer to dry it.
- Do not use a steam cleaner to clean it.

CAUTION (!):

- In case of damage or malfunction of your product DOT NOT use the product and please contact the authorized and trained representative.
- Do not use the product for a purpose not indicated in this manual
- MOBIAK S.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- MOBIAK S.A reserves the right to change the information contained in this document without previous notice
- There is a $\pm 3\%$. Deviation for the Dimensions of the Products.

GENERAL SAFETY WARNINGS:

- If you don't read this user manual, it is preferable not to use this product or another available part. If you don't understand the usage or the precautions, please contact with the dealer or the proper technical person before you use the product because may be caused damage.
- Keep away the product from heat sources.
- Save this document for future reference.
- Be careful when children are nearby and do not allow the children to play with the product.
- Do not attempt to lift the product by any removable parts.
- Do not use the product for a different purpose from the intended one.

WARRANTY:

The product is guaranteed for a period of one (1) year from the purchasing date.

The warranty apply on factory defects & DOES NOT apply on damages resulting from misuse, poor maintenance, modification, overuse, or non-compliance with the product's operating instructions.

Also DOES NOT apply on the parts, that can wear out during use or over time. Damages or defects caused by natural disasters, unauthorized maintenance or repair, power supply problems (where applicable) & carriers are not covered by the warranty.

No service or spare part is applied unless it has been inspected by an authorized service department or the manufacturer of the product. Shipping costs for goods of warranty are to be paid by the end user- customer or by the dealer.

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